



**Where to use :** With all sampling valves for a completely aseptic system.  
**How to use:** Manually-Instructions available  
**Use for:** 500ml Aseptic sampling bottle

**Accessories:**  
**Valve head:** W9 Q head 855543 (included)  
**Connection hose:** PTFE (included)  
**Filter:** Sterile (included)

**Material:**  
**Steel Parts:** AISI 304 / 316 / 316L(1.4404)  
**Membranes** FDA approved PTFE  
**O-rings:** Silicone and EPDM  
**Gasket:** PTFE and Silicone  
**Foot and stand:** AISI 316L and POM-C

**Key measurements**

<b>Drawing part A:</b>	250 mm.	9.84 In.
<b>Drawing part B:</b>	322 mm.	12.68 In.
<b>Bottle thread</b>	ISO GL45	
<b>Net weight:</b>	2.8 kg./pc	6.17 lb
<b>Temp. range:</b>	1 <sup>0</sup> C - 130 <sup>0</sup> C	34°F - 266°F
<b>Inlet pressure range:</b>	0 – 6 bar(g)	0 – 87 psi(g)
<b>Bottle pressure range:</b>	0-1,2 bar(g)	0-18 psi(g)