

Keofitt a/s  
Lumbyvej 11 H  
DK-5000 Odense C  
Denmark

Phone +45 63 16 70 80  
Fax +45 63 16 70 81  
E-mail: keofitt@keofitt.dk  
Web site: www.keofitt.dk

a/s reg. no.: 20 62 22  
VAT No. DK 42888214 (Import)

VAT No. DK 17112988 (Export)

## **WHY ON EARTH WOULD KEOFITT INTRODUCE AN INEXPENSIVE AND SIMPLE SAMPLING VALVE?**

Simply because most sampling cocks available in the market today are not sanitary.

Furthermore, most cocks do not fulfil the user's requirements for cost-efficiency, choice between manual and pneumatic operation, proven technology and documentation.

If you need quality sampling valves that are CIP-able, yet not steam-sterilisable (SIP) like a traditional Keofitt valve, our new SIMPLEX valve gives you value for money!

## **SIMPLEX – WHY SHOULD SAMPLING BE COMPLEX?**

Traditional steam-sterilisable Keofitt valves are dead leg free and fully drainable, making them completely safe. The Simplex sanitary sampling valve is based on the well-known technology of the existing W9 sterile sampling valve and is intended for sampling for physical and chemical analysis only.

## **SO WHAT MAKES THE SIMPLEX VALVE UNIQUE?**

The Simplex valve is CIP-able and comes with EHEDG certification to prove it! Furthermore, the Simplex valve shares most connection specifications and spare parts with the traditional Keofitt W9 valve.

### **SIMPLEX, THE NEW CIP-ABLE SAMPLING VALVE, OFFERS:**

- 316L machined steel body
- Manual or pneumatic operation
- W9 valve FDA-approved membranes
- Completely drainable
- EHEDG CIP certification
- 3.1 material certification
- Surface finish: Ra ≤ 0.8 µm