



Where to use :	With all sampling valves for a completely aseptic system.	
How to use:	Manually-Instructions available	
Use for:	500ml Aseptic sampling bottle	
Accessories:		
Valve head:	W9 Q head 855543 (included)	
Connection hose:	PTFE (included)	
Filter:	Sterile (included)	
Material:		
Steel Parts:	AISI 304 / 316 / 316L(1.4404)	
Membranes	FDA approved PTFE	
O-rings:	Silicone and EPDM	
Gasket:	PTFE and Silicone	
Key measurements		
Drawing part A:	164 mm.	6.46 In.
Drawing part B:	293 mm.	11.53 In.
Bottle thread	ISO GL 45	
Net weight:	1.7 kg./pc	3.75 lb
Temp. range:	1 ⁰ C - 130 ⁰ C	34°F - 266°F
Inlet pressure range:	0 – 6 bar(g)	0 – 87 psi(g)
Bottle pressure range:	0-1,2 bar(g)	0-18 psi(g)