

## **KEOFITT Plastic sleeve for Aseptic Sampling System**

ART. NO. 260225

## 

When handling glass bottles, there is always a risk that it may break. The new plastic sleeve offers convenience and also protects the bottle during handling and transport.

The plastic sleeve is a convenient and safe stand/holder for the KEOFITT Aseptic System with 1000 & 2000 ml bottles.

The top part of the plastic sleeve (black) is autoclaved together with the Aseptic System and bottle and after autoclavation the bottom part of the plastic sleeve (transparent) is screwed on to the top part.

It is important that the transparent part of the plastic sleeve is not autoclaved as the material is PVC which will get brittle after autoclavation.

The plastic sleeve is also available for 500 ml bottles, article no 260125.

## FEATURES

CERTIFICATION

The top part is attached to the bottle and Aseptic System while the bottom part is screwed on to the top part after autoclavation of the Aseptic System.



SOME PEOPLE WOULDN'T GAMBLE WITH THEIR SAMPLE™

		Aseptic Sampling System 260001	Aseptic Sampling System 260002	cem & Sampling Bottle Sampling bottle head 260030	Sampling bottle head for W15 860030	Sampling bottle head for W25 870030
ASEPTIC SYSTEM & SAMPLING BOTTLES	• • • • • • • • • • • • • • • • • • • •	including: Aseptic System Valve Pastic Sleeve for 500 ml bottle (260125) PTFE tube (Quick connection, 0,5 m, 550002) Filling tube for 500 ml bottles (260170)	<ul> <li>including:</li> <li>Aseptic System Valve</li> <li>Plastic Sleeve for 1000/2000 ml bottle (260225)</li> <li>PTFE tube (Quick connection, 0,5 m, 550002)</li> <li>Filling tube for 2000 ml bottles (260370)</li> </ul>	<ul> <li>PTFE tube (Quick connection, 0.5m, 550002)</li> <li>Plastic Sleeve for 500 ml bottle (260325)</li> <li>Filling tube for 0,5 litre bottles (260170)</li> </ul>	Mini tri-clamp. including: bottle (260325) bottle (260325)	Clamp 1" including: • Plastic Sleeve for 500 ml bottle (260325)